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
**INSTRUCTIONS FOR USE FOR 1K202, 1K173 AND 1K203-P6**
**I - USE**

These components are convenient for patients whose weight (carried load included) does not exceed 125Kg.

**II - PACK CONTENTS**

- 1 lamination connector (4 short prongs for 1K202, 1K203-P6 and 3 long prongs for 1K173),
- 1 spacer for lamination,
- This instruction sheet.

**III - SHAPING THE LAMINATION CONNECTOR**


It is highly recommended to adapt the shape of the connector prongs according to the positive.

**IV - RECOMMENDATIONS FOR LAMINATION**

- A** - Cover positive (1) with one or more layers of glass or carbon fabric (2) (according to the desired weight and strength).
- B** - For reinforcement purpose, install glass or carbon disc (3) where the lamination connector will be placed.
- C** - In order to reinforce the assembly, make braids (5) of glass or carbon fabrics pass through the holes of lamination connector (4).
- D** - Position lamination connector (4) onto positive (1) and check perpendicularity.
- E** - Install second disc (6) of glass or carbon fabric. The pyramid must necessarily be placed through this disc.
- F** - Clip spacer (7) on lamination connector (4).
- G** - Finish covering with one or more layers of glass or carbon fabric (8). Enclose the assembly in PVA bag (9) and inject resin. Before resin becomes hard, put string (10) around groove (11) located at the base of the lamination connector prongs. This will ensure a better strength.
- H** - When the resin is hard enough, cut off material excess and remove the spacer.

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