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
- .1 ensemble manette assemblée (1) comprenant :
 - 1 socle de manette (2)
 - 1 tirette (3)
 - 1 couvercle (4)
 - 1 insert (5)
 - 1 câble (6)
 - 1 vis pour matière plastique (7)
 - 1 vis Fhc M3 x 8 (8)
 - 1 vis Fhc M3 x 16 (9)
- .1 gaine (10)
- .2 pontets (11)
- .4 vis de pontet (12)
- .1 butée adhésive (13)
- .1 vis Fhc M3 x 20 (14)
- .1 vis Fhc M3 x 12 (15)
- .6 rivets tubulaires (16)
- . la présente notice

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

Fig 4

I PACK CONTENTS

- . 1 assembled handle (1) including :
 - 1 handle base (2)
 - 1 drawing device (3)
 - 1 cover (4)
 - 1 insert (5)
 - 1 cable (6)
 - 1 screw for plastics (7)
 - 1 screw M3 x 8 (8)
 - 1 screw M3 x 16 (9)
- . 1 sheath (10)
- . 2 cable retainers (11)
- . 4 screws for cable retainers (12)
- . 1 adhesive bearing (13)
- . 1 screw M3 x 20 (14)
- . 1 screw M3 x 12 (15)
- . 6 tubular rivets (16)
- . this instruction sheet

II PRODUCT CHARACTERISTICS

This handle is delivered in its version « automatically kept in unlocked position ». For an handle with automatic return, please see chapter IV.

Maximal release travel in version « automatically kept in unlocked position » : 14mm

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III WARNING


During some of the implementations described below, it may be necessary to disassemble the handle mechanism. In this case, the following instructions for assembling and disassembling the mechanism must be followed carefully :

- Disassemble the handle (fig.1) by unscrewing screw **8**, then **9** and then **7** (to reach screw **8**, the handle must be in low position, so that a 2mm hex wrench can be insert in hole **4a**).
- To reassemble the mechanism, proceed the same way by screwing screw **8**, then **7** and then **9**.

CAUTION ! : The drawing device 3 must be necessarily in high position each time you operate the handle. If it is in low position, screw 8 may move inside the mechanism and become unaccessible. In such a case, disassemble the mechanism, position screw 8 properly, then re-assemble the mechanism as indicated above.


IV HOW TO GET AN AUTOMATIC RETURN OF THE HANDLE

- a) Disassemble the handle as indicated in chapter III
- b) Raise flap **3a** with a screwdriver
- c) Level hooks **3b**
- d) Remove with a cutter the possible burrs that could hinder the mechanism operations
- e) Reassemble the mechanism.

IV FIXING THE HANDLE

WARNING ! This handle has been designed to be installed with two fixing screws. Any installation with just one screw could cause unvoluntary unlocking of the mechanism to which it is connected.

A) ON A LAMINATED SOCKET

- Make on the positive model a flat area not smaller than insert **5** (40mm height x 12 width)
- Position insert **5** in the middle of this area (tapped holes downwards as shown on the above diagram) by catching it under the two external fabrics layers.
- Laminate as usually.
- When the resin is completely hard, free the access to both tapped holes of insert **5**.
- Cut to length screws **8** and **9** (if necessary according to the socket wall thickness, replace them by screws **14** and **15** of the pack). **To avoid any damage of plastic parts, you must necessarily disassemble the mechanism** (chapter III) **before cutting the screws to length.**
- With a 2 mm hex wrench, assemble the handle by screwing screws **8** and **9** in both holes of insert **5** : recommended tightening torque 0.4 N.m (**tightening screw 9 too hard may jam the mechanism**)
Reminder : screwing screw **8** is possible only when the handle is in low position, so that a 2 mm hex wrench can be inserted in hole **4a**. Secure screws **8** and **9** with glue ref XC047.
- Determine the useful length of sheath **10** and cut it to length.
- Install sheath **10** on cable **6**
- Fasten sheath **10** with cable retainers **11** and screws **12** (or with rivets **16** according to the socket thickness).
- Then the cable can be connected to the knee. **WARNING ! Check that the release travel of the knee lock is not maximal when the handle is in high position.**
If the patient can hardly find the release button **4b** under his clothes, stick the adhesive bearing **13** on this button. We highly recommend to reinforce sticking of bearing **13** on button **4b** with PROTEOR XC050 glue.

B) ON A THERMOFORMED SOCKET


- Idem chapter **A**
- Position insert **5** in the middle of the flat area (tapped holes downwards as shown on the above diagram), taking care to put between the insert and the positive model a 2 mm thick wedge that has been cut to the outside dimensions of insert **5**, so that the insert be well embedded into the thermoformed material.
- Thermoform as usually
- After cooling, free the access to both tapped holes of the insert
- e), f), g), h), i), j), k) : Idem chapter **A**

V CABLE REPLACEMENT

In case of cable breakage, it can be replaced (cable ref : 3A24 25) as follows :

- The handle being positioned on the socket, fully unscrew screw **9** (do not act on screw **8**)
- Raise slightly the lower part of cover **4**
- Remove damaged cable **6**
- Place a new cable
- Screw screw **9**

Fig 4

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