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ADJUSTING PROCEDURES
FOR INTELLIGENT PROSTHESIS KNEE JOINT
IP PROGRAMMING UNIT : 1P301

I - KEY EXPLANATION


<RESET>
resets the micro-computer function and compensating for the 0 position of the needle valve

<SLOW>
slow walking speed

<DOWN>
decreases the data

<POWER>
turns on and off, and changes the mode

<CABLE>
2m (with two connectors on both ends)

<MID>
medium walking speed

<FAST>
fast walking speed

<UP>
increases the data

<WALKING SPEED>
measures walking speed and calculates average speed

<CLEAR>
clears the data and cancels a multistage setting

a. adjustment mode

Display	Operation procedure	SLOW	MID	FAST	UP/ DOWN	POWER	WALKING SPEED
	Connect the Adjuster to the knee joint, press the "POWER", then release it immediately.					○	
HELLO ! SP-51 V1.0	Wait for approximately three seconds.						
WAIT							
ADJUSTING MODE	Select the "MID" for normal walking.		○				
SELECT MV = 10	Perform a trial walk at middle (normal) speed, and select the optimum SV value for the knee swing using the "UP" and "DOWN" (Initial value : 10).				○		
STEPS = 7 MT = ---	Press the « Walking Speed » several times, and select the « STEPS » (5 to 7).						○
GO = 0 MT = 120	Walk a middle speed until "GO=0" is shown, and measure the speed MT. (7steps at prosthesis).						
V ? > 10 > ? T ? > 120 > ?	Check the value of middle speed, and select "SLOW" for slow walking.	○					
SELECT SV = 15	Perform a trial walk at slow speed, and select the optimum SV value for the knee swing using the "UP" and "DOWN".				○		
STEPS = 7 ST = ---	Press the "Walking Speed » several times, and select the "STEPS".						○
GO = 0 ST = 160	Walk at slow speed until "GO = 0" is shown, and measure the speed ST.						
V 15 > 10 > ? T160>120 > ?	Select « FAST » for fast walking.			○			
SELECT FV = 05	Perform a trial walk at fast speed, and select the optimum FV value for the knee swing using the "UP" and "DOWN".				○		
STEPS = 7 FT = ---	Press the "Walking Speed" several times and select the "STEPS".						○
GO = 0 FT = 80	Walk at fast speed until "GO = 0" is shown, and measure the speed ST,FT.						
V 15 > 10 > 05 T 160>120>80	When this message appears, press the "POWER" key					○	
GOOD-BYE NABCO	Multistage data is calculated on the basis of the détection data, and is transferred to the circuit board on the knee joint.						

b. Confirmation mode

Display	Operation procedure	SLOW	MID	FAST	UP/ DOWN	POWER	WALKING SPEED
	Connect the Adjuster to the knee joint, and press and hold the "POWER".						
HELLO ! SP-51 V1.0	Wait for approximately one second.					○	
CONFIRM MODE	When this message appears, release the "POWER".					○	
WAIT	Wait for approximately two seconds.					○	
CONFIRM MODE	When this message appears, perform a trial walk.					○	
CONFIRM MODE V5 = 15	The value on the second line varies with the walking speed.				○	○	
V1 = 05 B1 = 70	Press the "UP" and "DOWN" to display the multistage data registered to the knee joint.				○	○	
CONFIRM MODE V5 = 15	Keep the current status for three seconds. The previous message will reappear. When the "POWER" is pressed under status				○	○	
GOOD-BYE NABCO	the power turns OFF.					○	

b. Manual mode

Display	Operation procedure	SLOW	MID	FAST	UP/ DOWN	POWER	WALKING SPEED
	Connect the Adjuster to the knee joint, and press and hold the "POWER".						
HELLO ! SP-51 V1.0	Wait for approximately one second.					○	
CONFIRM MODE	While this message is being displayed, keep it pressed for another two seconds.					○	
MANUAL MODE	When this message appears, release the "POWER".					○	
WAIT	Wait for approximately two seconds.					○	
MANUAL MODE V = 15	The value on the second line shows the current fixed valve opening. press the "UP" and "DOWN" to select the optimum valve opening. When the "POWER" is pressed after the selection				○	○	
GOOD-BYE NABCO	with the valve opening of knee joint fixed, the power turns OFF.					○	

III - 1 : EROR MESSAGES

Message displayed : BATTERY VOLT. LOW

When : Turning on the power.

Cause : Exhausted battery.

Solution : Replace battery.

Special Note : this message will be displayed when voltage is low. This message may not appear if the device has been left on for a long period of time.

Message displayed : COM.ERROR !

When : When power is turned on during adjustment.

Cause : A) Exhausted battery.

B) Communication error between the programming unit and the circuit board.

1. Incorrect cable used.
2. Connection error.
3. Cable breakage.
4. Circuit board connected improperly.
5. Programming unit failure.
6. Dust or dirt on contact area.

Solution : A) Replace battery.

- B)
1. Use proper cable supplied.
 2. Insert connector fully.
 3. Replace cable.
 4. Press RESET key
 5. Replace programming unit.
 6. Clean contact area and keep plastic plugs in place.

Special Note : if error message disappears resume normal use.

If "COMM.ERROR" continues to be displayed and the power is shut off automatically, you need to press RESET once and turns the power on.

Message displayed : ST<MT ERROR !

When : After measuring ST data in adjusting mode.

Cause : The order of walking data is not SLOW, and is greater than or equal to MID.

Solution : Correct the order of SLOW walking modes.

Special Note : If the problem cannot be solved by the solution, the walking data of MID may be toolarge. Press CLEAR twice in succession to erase all the basic data. Then, load the basic data of MID again.

Message displayed : MT<FT ERROR !

When : After measuring ST data in adjusting mode.

Cause : The order of walking data is not MID, and is greater than or equal to FAST.

Solution : Correct the order of FAST walking modes.

Special Note : If the problem cannot be solved by the solution, the walking data of MID may be toolarge. Press CLEAR twice in succession to erase all the basic data. Then, load the basic data of MID again.

III - 2 : TROUBLESHOOTINGS

Use the following table if there is any trouble incurred during any phase of use.

Problem : **No display appears when power supply is turned on.**

Cause : A) Battery not connected to unit properly.
B) Battery exhausted.

Solution : A) Connect battery correctly.
B) Replace battery.

Problem : **A) The display blinks.**
B) Abnormal marks appear.
C) Display becomes weaker or disappears.

Cause : Exhausted battery.

Solution : Replace battery.

Problem : **A) Display blakens and is difficult to read.**
B) Display is faint and difficult to read.

Cause : A) Unit has been subjected to high or low temperatures for a long period of time.
B) There is an adjustment error of the unit.

Solution : A) Place in room with temperature between 10°C to 30°C for 1 hour.
Replace battery if unit does not recover in that time period.
B) Replace programming unit.

Problem : **A) Display appears only while the power key is being pressed, and disappears when released.**

Cause : A) Battery in the knee unit is not connected.
B) Battery in the knee unit is exhausted.

Solution : A) Connect battery in knee unit.
B) Replace battery unit in knee unit.

If display reads "BATTERY VOLT. LOW" or "COM. ERROR", refer to the section of Error Messages.

Problems during "MANUAL MODE" and/or "ADJUSTING MODE"

Problem : **Swing speed of the prosthesis does not change when data is changed.**

Cause : A) Battery is not connected.
B) "0" position error of needle valve.
C) Battery in knee unit is exhausted.
D) Cylinder module circuit board, motor, or needle valve.

Solution : A) Connect battery to motor.
B) Press the RESET key. Press the again key again if the first pressing not effective.
C) Replace battery in knee unit.
D) Return the knee unit to PROTEOR.

Problem : The remaining step indicator does not count down when walking speeds are measured.

- Cause : A) Proximity switch not connected.
 B) Damage to proximity switch, connection, or magnet.
- Solution : A) Return the knee unit to PROTEOR.
 B) Return the knee unit to PROTEOR.

If the display reads "ST<MT ERROR!" or "MT<FT ERROR!", refer to the section for ERROR MESSAGES.

Problem : Terminal Impact is too hard.

- Cause : Adjustment set too soft.
- Solution : Turn terminal impact adjustment screw clockwise to increase the dampening, however knee must still reach extension.

Problem : The prosthetic knee can reach full extension and therefore does not swing fully.

- Cause : Adjustment set too hard.
- Solution : Turn terminal impact adjustment screw counterclockwise to obtain full extension.

Problem : The knee unit makes grinding noise during swing phase.

- Cause : A) Rubber grommets are worn.
 B) Attachment bolts are loosened.
- Solution : A) Replace the rubber grommets.
 B) Retighten the fastening screws.

Problem : The swing of the prosthesis does not coordinate with a speed change.

- Cause : A) Connection to the battery, motor or proximity switch is loose.
 B) Battery in knee unit is exhausted.
 C) Damage to proximity switch or magnet.
 D) Pneumatic cylinder faulty.
- Solution : A) Insert connector fully.
 B) Replace battery in knee unit.
 C) Send the knee unit to PROTEOR for repair.
 D) Send the knee unit to PROTEOR.

Problem : Noise at any time.

- Cause : A) Wear or damage to rubber extension stops.
 B) Fastening bolts are loose.
 C) Wear of knee axis pin or bushings.
 D) Wear of cylinder pins or bushings.
 E) Wear of clevis bolt.
- Solution : A) Replace rubber extension stops.
 B) Retighten the loose bolts.
 C, D, E) Consult PROTEOR. ■