

# LAMES : 1A900

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
# BLADES : 1A900 CONNECTORS : 1K61, 1K62, 1K63 and 1K64

(Assembly and adjustment instructions for the CPOs)

## I FEATURES AND ADVANTAGES

1A900 : Large blade

- Mainly used for BK amputation or agenesia
- Flexible and dynamic
- Blade in composite material.

1K61 : Final connector that will be fixed under the socket

- Sturdy and easy to install
- Allows tilt adjustments in the frontal plane, as well as tilt and translation adjustments in the sagittal plane thanks to 1K60 or 1K66 pyramid connector (ordered separately).

1K62 : Temporary connector that will be fixed under the socket

- Sturdy and easy to install
- Allows tilt adjustments in the frontal plane, as well as tilt and translation adjustments in the sagittal plane thanks to 1K60 or 1K66 pyramid connector (ordered separately).


- Allows height adjustment of the prosthesis
- Fitted with 2 drill bushes to pre-drill the blade

1K63 : Final connector that will be laminated on the rear part of the socket

- Rigid, lightweight and very sturdy fixation


1K64 : Temporary connector that will be laminated on the rear part of the socket

- Allows tilt adjustments in the frontal plane, as well as height adjustment of the prosthesis
- Fitted with 2 drill bushes to pre-drill the blade

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## II GENERAL SAFETY INSTRUCTIONS

	<ul style="list-style-type: none"> <li>The product should be assembled only to the authorized prosthetic components mentioned in paragraph « V AUTHORIZED COMPONENTS ».</li> <li>The instructions for maximal weight and activity (running, footing) when using the product should be respected.</li> <li>Check that the top of the blade may not injure the patient and is not in contact with the socket.</li> <li>Do not heat the blade with a hot air gun, do not machine the blade to reduce its section, do not sand the sole too much (risk of removing the upper carbon layer).</li> <li>Do not drill the blade out of the fixing holes.</li> <li>Install the connectors only on the right part of the blade.</li> <li>Always use the clamping plates when fixing the connector to the blade.</li> <li>Do not forget to give to the patient the caring instructions sheet.</li> </ul>
---	---

## III CONTENT OF PACKAGING

<b>1A900 :</b> One blade One sole to be stucked under the blade One PROTEOR sticker	<b>1K61 :</b> One connector 4 screws FHC Ø6 2 screws H Ø10 One supporting plate	<b>1K62 :</b> One temporary connector 4 screws FHC Ø6 4 headless screws HC Ø6 One supporting plate A plate of abrasive grid to be cut.	<b>1K63 :</b> A lamination support 2 screws H Ø10 One supporting plate	<b>1K64 :</b> One temporary lamination support 2 screws CHC Ø8 2 washers One clamp A plate of abrasive grid to be cut.
--	---	---	---	---

## IV USE

These products are intended exclusively for lower limb prosthesis.

Blades 1A900 : Module selection depending on the patient weight and activity.							
Weight	40-49 kg (88-108 lb)	50-59 kg (110-130 lb)	60-69 kg (132-152 lb)	70-79 kg (154-174 lb)	80-89 kg (176-196 lb)	90-99 Kg (198-218 lb)	100-110 kg (220-242 lb)
Activity							
Footing		M1	M2	M3	M4	M5	M6
Running	M1	M2	M3	M4	M5	M6	


## V AUTHORIZED COMPONENTS

The blade can be assembled to 1K62 or 1K64 components for trial fitting and to 1K61 or 1K63 for final assembly.

1K61 and 1K62 connectors can be used with the following components : 1K66 in association with 1K203-P6, 1K202 and 1K173 ; or 1K60 in association with 1K172, 1K178, 1K207-P6, 1K203-P6 1K173.

**VI BLADE PREPARATION**

Sand the surfaces by hand before sticking the sole on the blade (XC090 Neoprene adhesive).

**VII ASSEMBLY**

Note : For efficiency purpose we recommend to fix the connector as low as possible on the right side of the blade (for example as shown here).

Knee joint center  
(15 mm below the joint line)

**Preparation :**

Trace the build line and the knee joint center line (patient standing with knee in complete extension).


**Sagittal plane :**

The build line should go through the knee joint center marked on the socket and 60 to 70 mm behind the point where the blade gets in contact with the ground.

The knee-ground assembly height should be from 30 to 50 mm higher than that of the patient (shoe included).

**Frontal plane :**

The build line should go through the center of the socket and through the center of the blade on all its length (the blade should be vertical with an external rotation of a few degrees).

**VIII USE OF 1K62 and 1K61 CONNECTORS :**

**1- Assembly with the 1K62 temporary connector**


Cut and stick (with Neoprene adhesive XC090) 2 pieces of abrasive grid (46x77 mm) on the supporting plate and on the contact face of the connector

Assemble connector 1K60 or 1K66 under the socket.

Assemble connector 1K62 on the blade.

Assemble the 2 parts of the prosthesis in the best possible compliance with the static alignment described in § VII, then finalize the assembly with the possible adjustments.

Tighten the assembly with the indicated torques and test the dynamic alignment of the prosthesis (see § « X DYNAMIC ALIGNMENT »). Re-adjust the blade position (height, tilting), then pre-drill it with Ø 4mm using the drill bushes of the temporary connector.

**2- Assembly with the 1K61 final connector**


Re-drill the blade with Ø10.5 (in the Ø4mm guide holes).

Replace the temporary connector with the final connector.

If necessary use 1K60 connector to finalize the alignment adjustments.

Re-tighten the screws up to the indicated torques and secure them with XC048 adhesive.

Cut the extra height of the blade if necessary, keeping at least 10 mm above the upper edge of the connector.


## IX USE OF 1K64 AND 1K63 POSTERIOR CONNECTORS

### 1– Assembly with 1K64 temporary connector

Once the blade position is determined, assemble connector 1K64 on the blade. Protect threads and drill passages with plastiline (on the sticking side).

Caution : If a sealing sheath is used, the upper edge of the connector should be at least 50 mm below the rear edge of the socket.

Fill the space between the connector and the socket with some thick and not brittle putty. Once the putty is hard, disassemble the blade and remove the clamp from the connector, then temporarily laminate the support on the socket. Cut 2 pieces of abrasive grid (46 x 77 mm) and stick them with XC090 neoprene adhesive on the two contact areas of the connector. Re-assemble the connector and the blade, then tighten up to the indicated torque. Test the dynamic alignment of the prosthesis (see § « X DYNAMIC ALIGNMENT ») and adjust the blade position. Once the position is well defined, drill the blade with Ø4mm using the drill bushes of the clamp.

### 2– Assembly with 1K63 final connector

Place the prosthesis on the alignment device, then mark the position of the blade and the socket. Cut the laminate around the temporary connector, then remove it from the socket. Re-drill the blade with Ø10.5mm (in the Ø4mm guide holes), then assemble final connector 1K63 on the blade. Protect the threads (on the sticking side). Position the blade and the socket on the alignment device again. Stick the final connector on the socket using some thick and not brittle putty. Remove the blade, then laminate strongly to fix the final connector on the socket. Re-assemble the blade, tighten up to the indicated torques and secure with XC048 adhesive. If necessary, cut the extra blade keeping at least 10 mm above the upper edge of the connector.

## X DYNAMIC ALIGNMENT

1– Jump from one foot to the other : the blade must stand flat on the ground (frontal view).

2– Hopping : The blade should propel upwards (neither forwards, nor backwards). If needed the blade is displaced backwards or forwards. If this has no effect, the tilt of the blade should be modified. In every case the patient should feel that his knee is free and easy to control (no hyper-extension due to the blade standing on the ground).

3– Running trial : The aim is to limit the muscular and ligament stresses and to decrease the energy used by the patient. For this purpose the height as well as the frontal and sagittal tilt are adjusted according to the patient's feeling. We look for a symetic stride at the patient's running speed.

- Feeling of sinking too far, of excessive softness : The blade module is probably too flexible. You can check this while running with small steps on an upward slope : if energy return is not sufficient, the blade is too flexible.

Action : Change the module. For a finer adjustment the blade can be displaced forwards compared to the socket.

- Feeling of shocks, discomfort : the blade module is probably too stiff.

Action : Change the module. For a finer adjustment the blade can be displaced backwards compared to the prosthesis.

- Feeling that the prosthesis swings outwards or inwards : the blade is probably too far inwards or outwards.

Action : It should be tilted in the valgum or varum direction.

- Feeling that the contact is divided into two stages (soft, then stiff) : the blade probably presses on the heel too early .

Action : backwards tilt should be increased in the sagittal plane.

- Feeling of falling before contact (the shoulder falls) : the prosthesis is probably too short.

Action : the prosthesis height should be increased.

- Feeling it is difficult to achieve the stride : the prosthesis is probably too long.

Action : the prosthesis height should be decreased.

**Please do not forget to give to the patient the caring instructions sheet (1A90098).**

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