## LAME DYNASPORT 1A900

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# 1A900 DYNASPORT BLADE

*(Caring instructions for the patient)*

**You have chosen the 1A900 sport blade DYNASPORT. Congratulations ! To get the best advantages from this product and to guarantee its functions and long life, it is important that you follow these instructions for use and care. Please do not hesitate to contact your CPO for any additional information you may need.**

- There are specific risks due to lower limb amputation (patient falling, trophic disorders of the residual limb caused by the prosthesis socket, ...). For these reasons, the patient must be trained to running by the CPO, or a physiotherapist, or in a disabled sports club or a rehabilitation center.
- Running on a treadmill can be considered under human supervision only and never when the patient is alone.

Please adhere to the following safety instructions

- Respect the maximal weight indicated and the recommended activities (running, footing) when using your product.
- Stop using the product in case of obvious damage.
- If you fall, check the blade and the fixations.
- Stop using the product if its behaviour becomes suspect (noise, change in flexibility) or if you have any doubt.
- Check the complete prosthesis after each use (screwed link, no noticeable looseness)
- Clean the prosthesis after each use.
- If needed contact your CPO and take the necessary actions (repair, replacement, checking by the manufacturer).
- Carefully handle the product.
- Do not use the product without the supplied sole, or if it is unsticked or damaged.
- Do not modify the product (cutting, ...).
- If you have to disassemble the blade, when re-assembling it adhere to the indicated tightening torques of the screws and stick them with PROTEOR XC048 adhesive).
- Do not place the prosthesis near a flame or an important source of heat (heater, ...).
- The rubber sole is a wear part. Please change it if it deteriorates to avoid damaging the blade (ref. 1A90002).

## 1A900 blades : Module selection depending on the patient weight and activity

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PROTEOR S.A.S.